

# **Dr. Maisch**

Any Column, Any Size, Any Media



## **ReproSil-XR**

eXtRa High Purity Phases

**MADE BY DR. MAISCH**

# CONTENT

- P 4 - eXtRa HIGH QUALITY HPLC PHASES
- P 5 - REPRODUCIBILITY AND PERFORMANCE
- P 6 - 9 - APPLICATIONS
- P 10 - 11 - REPOSIL XR C18MS



**REPOSIL-XR  
MADE BY DR. MAISCH**

From one of the biggest **High-Performance Liquid Chromatography (HPLC/UPLC)** - Column Manufacturers in Europe.

Choose Wisely - ReproSil-XR

Dr. Maisch HPLC GmbH sets a new standard in HPLC columns. We make everything - from the silica particle to the finished product. Our streamlined manufacturing process allows delivery of the highest performance at an exceptional value!


eXtRa high quality HPLC phases

Extra High purity silica and uniform bonded phase coverage translates to symmetrical peaks for acids/bases, and predictable reversed phase selectivity. Whether routine analysis or new method developments, use ReproSil XR columns to get premium performance for a low budget.


ReproSil-XR	Specifications
Modifications	Silica, C18, C18(MS), C8, Phenyl, CN, SCX, SAX
Endcapping	(for C18, C8, Phenyl)
Surface area	250 m <sup>2</sup> /g , 100 m <sup>2</sup> /g
Pore size	120 Å, 300 Å
Particle size	1.5, 3, 3.5, 5, 7 and 10 µm

Reproducible Methods Start with Reproducible Columns

Highly Reproducible Selectivity




Tightly controlled silica synthesis and carefully performed bonding result in minimum variations of selectivity and capacity factors.




Consistent batch-to-batch capacity for acidic, basic & neutral compounds.

Expect competitive performance




Efficient manufacturing processes are the driving force for competitive pricing. But certainly not at the risk of competitive performance. ReproSil-XR columns show similar or better efficiency and asymmetry for challenging basic and acidic compounds compared to industry leading LC columns.



ReproSil-XR columns compares favorable to leading competitor columns.

Amoxicillin (U.S.P.)



**Column:** ReproSil-XR 120 C18, 5 µm, 120 Å, 250 x 4.6 mm, (PN: rx15.9e.s2546)


**Mobile phase:** Isocratic, 50 mM KH<sub>2</sub>PO<sub>4</sub>, pH 5 : MeCN 96:4 (v:v)

**Flow rate:** 1.5 ml/min

**Detector:** UV at 230 nm


**Column temp:** 30° C

**Injection vol:** 10 µl



	ReproSil-XR C18	Comp P C18	Comp G C18	Comp H C18	Comp K C18
Theoretical plates	8650	7240	6230	7060	7680
Tailing factor	1.0	1.2	1.1	1.0	1.0

Standards Mix




- 1. Nicotine
- 2. Quinine
- 3. Diphenhydramine
- 4. Nortriptyline
- 5. 2,4 -Dichlorophenoxy acetic acid
- 6. 2,4 -Dichlorophenoxy propionic acid
- 7. 2,4 -Dichlorophenoxy butyric acid

**Column:** **ReproSil-XR** 120 C18, 5 µm,  
120 Å,  
150 x 4.6 mm  
(PN: rx15.9e.s1546)

**Mobile phase:**  
MeCN : 50 mM KH<sub>2</sub>PO<sub>4</sub>  
50:50 (v:v) pH 3 at 35°C

Baseline resolution challenging separations.

Organic Acids




- 1. 5.5 mg/ml Tartaric acid
- 2. 11.8 mg/ml Malic acid
- 3. 15.2 mg/ml Lactic acid
- 4. 0.13 mg/ml Maleic acid
- 5. 0.18 mg/ml Fumaric acid

**Column:** **ReproSil-XR** 120 C18, 5 µm,  
120 Å,  
150 x 4.6 mm  
(PN: rx15.9e.s1546)

**Mobile phase:** 25 mM KH<sub>2</sub>PO<sub>4</sub> : MeOH  
97:3 (v:v)

**Flow rate:** 0.7 ml/min  
**Detector:** UV at 220 nm  
**Injection vol:** 1 µl

Acids, Bases and Neutrals




- 1. Uracil
- 2. Amitriptyline
- 3. Phenylvaleric acid
- 4. Toluene

**Column:** **ReproSil-XR** 120 C18, 5 µm,  
120 Å,  
150 x 4.6 mm  
(PN: rx15.9e.s1546)

**Mobile phase:**  
MeCN : 50 mM KH<sub>2</sub>PO<sub>4</sub>  
50:50 (v:v) pH 3 at 35°C

Symmetric peaks for acids and bases.

Carbamate & Urea




- 1. Fenuron
- 2. Monuron
- 3. Fluometuron
- 4. Diuron
- 5. Methiocarb
- 6. Linuron
- 7. Neburon

**Column:** **ReproSil-XR** 120 C18, 5 µm,  
120 Å,  
150 x 4.6 mm  
(PN: rx15.9e.s1546)

**Mobile phase:** **A:** 25 mM KH<sub>2</sub>PO<sub>4</sub>  
pH 3.2  
**B:** MeCN

**Gradient:** (Time, %B): (0,30), (20,60)  
**Flow rate:** 1.0 ml/min  
**Detector:** UV at 240 nm  
**Injection vol:** 10 µl

Water Soluble Vitamins



**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 150 x 4.6 mm (PN: rx15.9e.s1546)

**Mobile phase:** **A:** 100 mM NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub> Buffer, pH 4.5  
**B:** MeCN


**Gradient:** (Time, %B): (0,5), (4,5), (10,40), (15,40)

**Flow rate:** 1.0 mL/min

**Detector:** UV at 254 nm

**Injection vol:** 2 µl

Folium Ginkgo




**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 150 x 4.6 mm (PN: rx15.9e.s1546)

**Mobile phase:** MeOH:THF:H<sub>2</sub>O 25:10:65 (v:v:v)

**Flow rate:** 1.0 ml/min

**Detector:** ELSD

Phenols



**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 150 x 4.6 mm (PN: rx15.9e.s1546)

**Mobile phase:** **A:** 1% CH<sub>3</sub>COOH in H<sub>2</sub>O  
**B:** 1% CH<sub>3</sub>COOH in MeOH


**Gradient:** (Time, %B): (0,45), (6,45), (8,60), (15,80)

**Flow rate:** 1.0 mL/min

**Detector:** UV at 280 nm

**Injection vol:** 10 µl

Simvastatin Tablets (U.S.P)



**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 250 x 4.6 mm (PN: rx15.9e.s2546)

**Mobile phase:** 0.39% KH<sub>2</sub>PO<sub>4</sub> pH 4.5 : MeCN 35:65 (v:v)


**Flow rate:** 1.5 ml/min

**Detector:** UV at 238 nm

**Column temp:** 45°C

**Injection vol:** 10 µl

Paracetamol effervescent tablets on chinese pharmacopeia assay




**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 250 x 4.6 mm (PN: rx15.9e.s2546)

**Mobile phase:** KH<sub>2</sub>PO<sub>4</sub> : MeOH 80:20 (v:v)

**Flow rate:** 1.0 mL/min

**Detector:** UV at 254 nm

Ibuprofen (U.S.P.)



**Column:** **ReproSil-XR** 120 C18, 5 µm, 120 Å, 150 x 4.6 mm (PN: rx15.9e.s1546)

**Mobile phase:** H<sub>2</sub>O (adjusted with H<sub>3</sub>PO<sub>4</sub> to a pH of 2.5) : MeCN 65:35 (v:v)

**Flow rate:** 2 ml/min

**Detector:** UV at 214 nm

**Injection vol:** 5 µl




## Difficult acids/bases or chelates? No Problem!

### Outstanding features of ReproSil-XR

<b>High Polarity:</b>	Different selectivity to current phases.
<b>High Base Retention:</b>	Uses high % organic: helps MS detection.
<b>Acids, Bases, Chelates:</b>	All the benefits of pure silica.
<b>100% Aqueous Stability:</b>	Resists phase collapse.
<b>Compatible with MS buffers:</b>	Good with volatile buffers.
<b>High Stability:</b>	No embedded phase to limit stability.
<b>Low MS bleed:</b>	Helps MS Detection.


### Separate Highly Basic Components at Neutral PH



1. Phenylephrine
2. Diphenhydramine
3. Amitriptyline

**Column:** ReproSil-XR 120 C18MS, 1.5 μm, 120 Å, 50 x 2.0 mm  
**Mobile phase:** (PN: rx115.9ms.s0502) 50mM NH<sub>4</sub>COOH, pH 7: MeOH (20:80)  
**Column Temp:** 40°C  
**Flow rate:** 0.2mL/min  
**Detector:** UV at 210nm

### Excellent Peak Shape for Acids, Bases and Chelates




- |                              |                    |                       |
|------------------------------|--------------------|-----------------------|
| 1. 2,4-Dichlorophenoxyacetic | 1. PsuedoEphidrine | 1. Pyridine,          |
| 2. Propionic acid            | 2. Diphenhydramine | 2. 8-Hydroxyquinoline |
| 3. Butyric acids             | 3. Amitriptyline   | 3. 2,2' Dipyridyl     |

**Column:** ReproSil-XR 120 C18MS, 5 μm, 120 Å, 150 x 4.6 mm  
**Mobile phase:** (PN: rx15.9ms.s1546) 50% MeCN : 50% 50 mM  
**Column Temp:** KH<sub>2</sub>PO<sub>4</sub> pH 3,  
**Flow rate:** 35°C

## MS Bleeding reduced by Factor 17

### Comparison of MS Bleeding




ReproSil-XR 120 C18MS shows minimal bleed when used in LC/MS with TFA gradients. The first trace shows bleed levels typical of polar embedded/encapped product.

### 4 Important Benefits

- Higher methanol level for improved MS sensitivity
- 2.5 x lower ammonium formate entering the MS (10% aqueous vs. 25%)
- Half the backpressure due to lower viscosity.
- Faster analysis due to the potential to use higher flowrate without excessive backpressure.

### Pesticides



Fragmentor Voltage 150 Collision Energy 45 Ionization Mode ESI

**Column:** ReproSil-XR 120 C18MS, 1.5 μm, 120 Å, 10 x 2mm  
**Mobile phase:** A: H<sub>2</sub>O  
 B: MeCN  
**Gradient:**

A:	0	2	3	5	6	8	8.1	10
B:	5	5	35	35	95	95	5	5

**Flow rate:** 0.4ml/min  
**Detector:** 6400 Series Triple Quadrupole