

Inertsil® AX


Physical Properties

- Silica : 3 Series High Purity Silica Gel
- Particle Size : 5 μm
- Surface Area : 450 m^2/g
- Pore Size : 100 \AA (10 nm)
- Pore Volume : 1.05 mL/g
- Bonded Phase : Diethylaminopropyl Groups
- End-capping : None
- Carbon Loading : 17 %
- AEC : 0.4 meq/g
- USP Code : -
- pH Range : 2 ~ 7.5




Inertsil AX has diethylamino groups bonded to silica gel by an alkyl chain. The diethylamino groups offer anionic functions required for anion exchange chromatography. It is mainly used for analyses of acidic compounds. Conventional ion-exchange columns used to show inconsistent results from lot to lot. However, Inertsil AX is manufactured under strict quality control in order to offer excellent lot-to-lot reproducibility. The retentivity of Inertsil AX is influenced by the concentration of buffer. The retention time can be adjusted by the concentration of buffer (Refer to Fig. 2).

Figure 1 : Comparison of lot-to-lot reproducibility with other brands




Conditions

Column Size : 5 μm , 150 \times 4.6 mm I.D.
 Eluent : 60 mM KH_2PO_4 (pH 3.0, H_3PO_4)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 $^\circ\text{C}$
 Detection : UV 254 nm
 Injection Vol. : 1 μL
 Sample : 1. Cytidine 5'-monophosphate (CMP)
 2. Uridine 5'-monophosphate (UMP)




Cytidine 5'-monophosphate (CMP)



Uridine 5'-monophosphate (UMP)

Figure 2 : Effect of buffer concentration in eluent



Conditions

Column : Inertsil AX (5 μm , 150 \times 4.6 mm I.D.)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 $^\circ\text{C}$
 Detection : UV 254 nm
 Injection Vol. : 10 μL
 Sample : 1. Cytidine 5'-monophosphate (CMP)
 2. Adenine 5'-monophosphate (AMP)
 3. Uridine 5'-monophosphate (UMP)
 4. Guanosine 5'-monophosphate (GMP)
 5. Xanthosine 5'-monophosphate (XMP)

Analytical Columns

Particle Size: 5 µm	Length \ I.D. (mm)		1.0		1.5					
	33		5020-80111	5020-80121						
	50		5020-80112	5020-80122						
	75		5020-80113	5020-80123						
	100		5020-80114	5020-80124						
	150		5020-80115	5020-80125						
	250		5020-80116	5020-80126						
	Length \ I.D. (mm)		2.1		3.0		4.0		4.6	
	33		5020-07211	5020-07221	5020-07231	5020-07241				
	50		5020-07212	5020-07222	5020-07232	5020-07242				
75		5020-07213	5020-07223	5020-07233	5020-07243					
100		5020-07214	5020-07224	5020-07234	5020-07244					
150		5020-07215	5020-07225	5020-07235	5020-07245					
250		5020-07216	5020-07226	5020-07236	5020-07246					

* End-fittings are 1/16" Waters-compatible.

* For maximum operating pressure information, please refer to page 46.

Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)		Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)	
			Particle Size		Particle Size	
			5 µm		5 µm	
1.0	10	1.0	5020-19233		5020-19283	
1.5, 2.1		1.5	5020-19333		5020-19383	
2.1, 3.0		3.0	5020-19133		5020-19183	
4.0, 4.6		4.0	5020-19033		5020-19083	
2.1, 3.0	20	3.0	5020-19533		5020-19583	
4.0, 4.6		4.0	5020-19433		5020-19483	
Holder for Cartridge Guard Column E			For 10 mm Length		5020-08500	
			For 20 mm Length		5020-08550	

* End-fittings are 1/16" Waters-compatible.

* For maximum operating pressure information, please refer to page 46.

Reversed Phase Columns

HILIC Columns

Normal Phase Columns

SEC Columns

Ion Exchange Columns

Application Specific Columns

Guard Columns

Preparative Columns

Capillary Columns

Applications

Cat. No. Index