


Fullerenes

"Bucky Balls"




CHROM 7210

1. C70
2. C60


Column: Jordi™ DVB, 500Å, 500 x 10mm (6) (Part No. 100567)
Mobile Phase: Trichlorobenzene
Flow Rate: 1.5mL/min
Column Temp: 145°C
Detector: RI

Buckminsterfullerenes



Column: Vydac® C4, 5µm, 250 x 10mm (Part No. 201TP510)
Mobile Phase: Acetonitrile:Toluene (50:50)
Flow Rate: 5mL/min

Regioisomeric Fullerene Derivatives




1. *trans*-1-C62(COOEt)₄
2. *trans*-2-C62(COOEt)₄
3. *trans*-3-C62(COOEt)₄
4. *trans*-4-C62(COOEt)₄
5. *e*-C62(COOEt)₄

Column: Grom™ Sil 100 Norm Ph-1 ST, 5µm, 250 x 4mm
Mobile Phase: Toluene
Flow Rate: 1.5mL/min
Detector: UV at 340nm

Cooking Oil

Used Cooking Oil




CHROM 7215

1. Polymerized
2. Polymerized
3. Triglycerides
4. MW 1800-3000
5. Triglycerides MW 930
6. Diglycerides MW 620
7. Monoglycerides MW 400
8. Free Fatty Acids MW 300

Column: Jordi™ DVB, 10³Å, 500Å, 500mm x 10mm (2) (Part No.s 100571, 100567)
Mobile Phase: Tetrahydrofuran
Flow Rate: 2.0mL/min
Detector: UV at 254nm

Surfactants

Microbore Analysis of Glucopyranosides



CHROM 10310

1. *n*-octyl Glucopyranoside (50ng on column)
2. *n*-decyl Glucopyranoside (50ng on column)
3. *n*-dodecyl Glucopyranoside (50ng on column)


Column: Alltech® Alltima™ C18, 5µm, 150 x 1mm (Part No. 88384)
Mobile Phase: Methanol:Water (90:10)
Flow Rate: 0.05mL/min
Detector: ELSD



Polymers

Polyurethane Curatives

CHROM 7472




- 1. Polytetramethylene Glycol
- 2. Ethechure 300

Column: Jordi™ DVB, 500Å, 500mm x 10mm (6)
(Part No. 100567)
Mobile Phase: Chloroform
Flow Rate: 1.5mL/min
Column Temp: 50°C
Detector: RI

Polystyrene

CHROM 7213




- 1. 8.42M
- 2. 1.3M
- 3. 791K
- 4. 354K
- 5. 98.9K
- 6. 43K
- 7. 10.1K
- 8. 2.6K
- 9. 530K
- 10. 120K

Column: Jordi™ DVB, Mixed Bed Linear, 500mm x 10mm
(Part No. 100583)
Mobile Phase: Tetrahydrofuran
Flow Rate: 1.5mL/min
Column Temp: 40°C
Detector: RI

miscellaneous

Polyethylene and Hydrocarbon Standards

CHROM 7220




- | | | |
|----------------|---------|---------|
| 1. PE-32000 MW | 14. C36 | 27. C10 |
| 2. PE-2155 MW | 15. C34 | 28. C9 |
| 3. C58 | 16. C32 | 29. C8 |
| 4. C56 | 17. C30 | 30. C7 |
| 5. C54 | 18. C28 | 31. C6 |
| 6. C52 | 19. C26 | 32. C5 |
| 7. C50 | 20. C24 | |
| 8. C48 | 21. C22 | |
| 9. C46 | 22. C20 | |
| 10. C44 | 23. C18 | |
| 11. C42 | 24. C16 | |
| 12. C40 | 25. C14 | |
| 13. C38 | 26. C12 | |

Column: Jordi™ DVB, 500Å, 500mm x 10mm (6)
(Part No. 100567)
Mobile Phase: Trichlorobenzene
Flow Rate: 1.5mL/min
Column Temp: 145°C
Detector: RI

Tween 80


CHROM 7809



Column: Alltech® Alltima™ C18, 5µm, 250 x 4.6mm
(Part No. 88056)
Mobile Phase: Water:Methanol (40:60)
Flow Rate: 1.0mL/min
Detector: ELSD

Dimethicone

CHROM 8636



Column: Alltech® Alltima™ C8, 5µm, 150 x 4.6mm (Part No. 88072)
Mobile Phase: A: Acetonitrile B: Chloroform
Gradient:

Time:	0	20
%B:	25	75

Flow Rate: 1.0mL/min
Detector: ELSD


more applications

View our complete searchable chromatogram database at www.discoverysciences.com/chromdb/

Polymers

Polyethylene Glycol (MW=600)

CHROM
8442



Column: Alltech® Adsorbosphere™ C8, 5µm, 250 x 4.6mm (Part No. 287112)
Mobile Phase: A: Methanol B: Water
Gradient:


Time:	0	45
%B:	60	20

Flow Rate: 0.6mL/min
Detector: ELSD

Phenones/Fast Gradient Column Formats

Phenone Homologs

CHROM
9467



1. Uracil
2. Theophylline
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Heptanophenone
9. Octanophenone


Columns: Alltech® Alltima™ C18, 3µm, 20 x 4.6mm Expedite™ (Part No. 43451)
Mobile Phase: A: 0.1% Trifluoroacetic Acid in Water
 B: 0.08% Trifluoroacetic Acid in Acetonitrile
Gradient:

Time:	0	2
%B:	5	100

Flow Rate: 2.0mL/min
Detector: UV at 254nm

Phenone Homologs

CHROM
9471



1. Uracil
2. Theophylline
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Heptanophenone
9. Octanophenone


Columns: Alltech® Alltima™ C18, 3µm, 10 x 2.1mm Expedite™ (Part No. 43836)
Mobile Phase: A: 0.1% Trifluoroacetic Acid in Water
 B: 0.08% Trifluoroacetic Acid in Acetonitrile
Gradient:

Time:	0	2
%B:	5	100

Flow Rate: 1.0mL/min
Detector: UV at 254nm

Phenone Homologs

CHROM
9468



1. Uracil
2. Theophylline
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Heptanophenone
9. Octanophenone


Columns: Alltech® Alltima™ C18, 3µm, 10 x 4.6mm Expedite™ (Part No. 43837)
Mobile Phase: A: 0.1% Trifluoroacetic Acid in Water
 B: 0.08% Trifluoroacetic Acid in Acetonitrile
Gradient:

Time:	0	2
%B:	5	100

Flow Rate: 2.0mL/min
Detector: UV at 254nm

Phenone Homologs

CHROM
9460



1. Uracil
2. Theophylline
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Heptanophenone
9. Octanophenone


Column: Alltech® Prevail™ C18, 3µm, 20 x 2.1mm, Expedite™ (Part No. 43827)
Mobile Phase: A: 0.1% Trifluoroacetic Acid in Water
 B: 0.08% Trifluoroacetic Acid in Acetonitrile
Gradient:

Time:	0	2	3
%B:	0	100	100

Flow Rate: 1.0mL/min
Detector: UV at 254nm

Phenone Homologs

CHROM
9462



1. Uracil
2. Theophylline
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Heptanophenone
9. Octanophenone

Column: Alltech® Alltima™ C18, 3µm, 20 x 2.1mm Expedite™ (Part No. 43803)
Mobile Phase: A: 0.1% Trifluoroacetic Acid in Water
 B: 0.08% Trifluoroacetic Acid in Acetonitrile
Gradient:


Time:	0	2
%B:	5	100

Flow Rate: 1.2mL/min
Detector: UV at 254nm

Biofuels

Triacylglycerols in Canola Oil

CHROM 10749



- 1. LLnLn
- 2. LLLn
- 3. OLnL
- 4. LLO
- 5. LOO
- 6. PLO
- 7. OOO
- 8. POO
- 9. POP
- 10. PPP


Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
 Mobile Phase: A: Acetonitrile B: Dichloromethane
 Gradient:

Time:	0	20	25	30	35
%B:	30	70	70	30	30

 Flow Rate: 1.0mL/min
 Detector: ELSD

Triacylglycerols in Red Palm Oil

CHROM 10750



- 1. PLL
- 2. LOO
- 3. PLO
- 4. PLP
- 5. MOP
- 6. OOO
- 7. POO
- 8. POP
- 9. SOO
- 10. PPP
- 11. POS
- 12. PPS
- 13. SOS


Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
 Mobile Phase: A: Acetonitrile B: Dichloromethane
 Gradient:

Time:	0	20	25	30	35
%B:	30	70	70	30	30

 Flow Rate: 1.0mL/min
 Detector: ELSD

Biodiesel in Spec with EN14105 Method

CHROM 10644



- 1. Mono-acylglycerol (MAG)
- 2. MAG
- 3. MAG
- 4. Fatty Acid Methyl Ester (FAME)
- 5. MAG
- 6. FAME
- 7. FAME
- 8. MAG
- 9. FAME
- 10. Di-acylglycerol (DAG)
- 11. Tri-acylglycerol (TAG)


Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
 Mobile Phase: A: Acetonitrile B: Dichloromethane
 Gradient:

Time:	0	5	30	32	35
%B:	0	15	70	70	0

 Flow Rate: 1.0mL/min
 Detector: ELSD

Triacylglycerols in Used Chinese Frying Oil

CHROM 10751



- 1. LLnLn
- 2. LLLn
- 3. LLL
- 4. OLLn
- 5. PLLn
- 6. LLO
- 7. PLL
- 8. LOO
- 9. LLS
- 10. PLO
- 11. PLP
- 12. OOO
- 13. SLO
- 14. POO
- 15. POP


Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
 Mobile Phase: A: Acetonitrile B: Dichloromethane
 Gradient:

Time:	0	20	25	30	35
%B:	30	70	70	30	30

 Flow Rate: 1.0mL/min
 Detector: ELSD

Triacylglycerols in Soy Oil

CHROM 10748



- 1. LLnLn
- 2. LLLn
- 3. LLL
- 4. OLLn
- 5. PLLn
- 6. LLO
- 7. PLL
- 8. LLO
- 9. LLS
- 10. PLO
- 11. PLP
- 12. OOO
- 13. SLO
- 14. POO
- 15. POP

Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
 Mobile Phase: A: Acetonitrile B: Dichloromethane
 Gradient:

Time:	0	20	25	30	35
%B:	30	70	70	30	30

 Flow Rate: 1.0mL/min
 Detector: ELSD


more applications

View our complete searchable chromatogram database at www.discoverysciences.com/chromdb/

Biofuels

Triacylglycerols in Lard

CHROM 10752



1. PPO
2. PPP
3. POS
4. SSS

Column: Alltech® Alltima™ HP C18 HiLoad, 5µm, 250 x 4.6mm (Part No. 87698)
Mobile Phase: A: Acetonitrile B: Dichloromethane
Gradient:


Time:	0	20	25	30	35
%B:	30	70	70	30	30

Flow Rate: 1.0mL/min
Detector: ELSD

Fire Retardants

Brominated Flame Retardants

CHROM 10261




1. Hexabromo Cyclododecane
2. Tetrabromo Bisphenol A (Allyl Ether)

Column: Alltech® Alltima™ C18, 5µm, 250 x 4.6mm (Part No. 88056)
Mobile Phase: Water:Acetonitrile (5:95)
Flow Rate: 1.0mL/min
Detector: ELSD

Standards

Standards Mix

CHROM 10696




1. Nicotine
2. Quinine
3. Diphenhydramine
4. Nortriptyline
5. 2,4, Dichlorophenoxy Acetic Acid
6. 2,4, Dichlorophenoxy Propionic Acid
7. 2,4, Dichlorophenoxy Butyric Acid

Column: GraceSmart™ C18, 5µm, 150 x 4.6mm (Part No. 5138812)
Mobile Phase: 50mM Potassium Phosphate, pH3.0:Acetonitrile (50:50)
Detector: UV at 254nm

Acids, Bases and Neutrals

CHROM 10695




1. Uracil
2. Amitriptyline
3. Phenylvaleric Acid
4. Toluene

Column: GraceSmart™ C18, 5µm, 150 x 4.6mm (Part No. 5138812)
Mobile Phase: 50mM Potassium Phosphate, pH3.0:Acetonitrile (50:50)
Detector: UV at 254nm

Dyes

Color Additives - Dyes



1. Tartrazine
2. Amaranth
3. Indigo Carmine
4. New Coccine
5. Sunset Yellow
6. Fast Green
7. Brilliant Blue
8. Erythosine
9. Acid Red
10. Phloxine
11. Rosé Bengale


Column: Grom™ Sil 120 ODS-5 ST, 5µm, 150 x 4mm (Part No. 5113642)
Mobile Phase: A: 10mM Ammonium Phosphate, pH6.0:Methanol (90:10)
 B: 10mM Ammonium Phosphate, pH6.0:Methanol (20:30)
Gradient:

Time:	0	40	60
%B:	0	100	100

Flow Rate: 0.8mL/min
Detector: UV at 254nm

Solvents

Caprolactams




1. N-Methylol Caprolactam
2. Caprolactam
3. N,N-Methylene-bis(Caprolactam)
4. N,N-oxy-bis(Methylene)-Caprolactam

Column: Jones Genesis®, 4µm, 150 x 4.6mm (Part No. FM15960E)
Mobile Phase: Acetonitrile:20mM Phosphate, pH6.8 (22:78)
Flow Rate: 1.0mL/min
Detector: UV at 205nm

Fuels

Petrochemicals




- 1. Aliphatics
- 2. Mono-aromatics
- 3. Di-aromatics
- 4. Polyaromatics

Column: Jones Genesis® Petro XL, 3µm, 250 x 4.6mm (Part No. FM25991E)
Mobile Phase: Hexane
Flow Rate: 1.0mL/min
Detector: UV at 210nm

Chelators

Metal Chelators

CHROM 8584




- 1. 2,2'-Dipyridyl
- 2. 8-Hydroxyquinoline

Column: Alltech® Alltima™ C18, 5µm, 150 x 4.6mm (Part No. 88052)
Mobile Phase: 0.025M K₂HPO₄, pH7.0:Tetrahydrofuran: Methanol (50:25:25)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

Paraffins

Paraffin Analysis

CHROM 10429




Column: Alltech® Alltima™ HP C18 HiLoad, 3µm, 100 x 4.6mm (Part No. 87685)
Mobile Phase: Acetonitrile:Dichloromethane (60:40)
Flow Rate: 1.0mL/min
Detector: ELSD

Cyclodextrins

Cyclodextrins

- 1. γ-Cyclodextrin
- 2. α-Cyclodextrin
- 3. β-Cyclodextrin



Column: Grom™ Sil 120 ODS-5 ST, 5µm, 150 x 4mm (Part No. 5113642)
Mobile Phase: Water:Methanol (95:5)
Flow Rate: 0.8mL/min
Detector: RI

miscellaneous