

Reveleris® iES System Purifies Ferulic Acid in Rhizoma Chuanxiong Extract with Greater Recovery

Reveleris® iES Flash Chromatography System

Introduction

Medicinal plants are an important resource for drug development. They are commonly analyzed in biochemical, pharmaceutical, and clinical research for their complex constituents and low quantities of bioactive ingredients.

Ferulic acid is known to be one of the main bioactive components found in Rhizoma chuanxiong.¹ The plant has been extracted here and used for the separation and identification of its scarce ingredients using the Reveleris® iES flash chromatography system.

Experimental

Extraction Conditions

Extract type: Pulverized

Weight: 2 g

Extraction solvent: Ether

Solvent volume: 40 mL

Ultra-sonication: 20 min

Run Conditions

Cartridge: Reveleris® silica 12g

Flow rate: 25 mL/min

Equilibration: 3.0 min

UV1 wavelength: 254 nm

UV2 wavelength: 280 nm

ELSD carrier: Iso-propanol

Injection type: Liquid

Solvent A: Hexane

Solvent B: Ethyl acetate

Gradient Method

Step	Time (min.)	%B
1	0	0
2	3	0
3	8.8	5
4	0.5	6
5	2.6	33
6	5.2	98
7	2.6	98

HPLC Conditions

Column: Alltima™ silica 5µm, 250 x 4.6 mm (PN: 88171)

Flow rate: 1 mL/min

UV Detector: 280 nm

Mobile phase: A: 1% Formic acid in Hexane

B: 1% Formic acid in Ethyl acetate

Gradient Method

Step	Time (min.)	%B
1	0	0
2	10	50
3	10	100



Results and Discussion




Figure 1: Separation of Rhizoma chuanxiong extract shows the detection of peaks using the UV and the evaporative light scattering detectors.




Figure 2: Ferulic acid recovery = 97% (calculated by peak area)

Conclusion

Using the sensitive Reveleris® iES flash chromatography system, the critical components of the plant extract at low levels have been isolated and collected. The RevealX™ technology also allows non-chromophoric peaks to be detected with the chromophoric ones in the same run

for fraction collection and further analysis. This allows purification of compounds present in the natural product extract, compounds that are not always easily detectable using traditional flash chromatography system and are present at low amounts and useful for medicinal purposes.

References

1. Chen, X.; Kong, L.; Su, X.; Fu, H.; Ni, J.; Zhao, R.; Zou, H.; Separation and identification of compounds in Rhizoma chuanxiong by comprehensive two-dimensional liquid chromatography coupled to mass spectrometry; J. Chromatography A, 1040, (2004), pp. 169-178.

Americas
Deerfield, Illinois
Tel: +1 847.948.8600
discoverysciences@
grace.com

Europe
Lokeren, Belgium
Tel: +32 (0)9.340.65.65
discoverysciences.EU@
grace.com

India
Pune
Tel: +91 20.6644.9900
pune@grace.com

China
Shanghai
Tel: +86 21.5467.4678
dsbiz.asia@grace.com

Australia
Rowville
Tel: +61 3.9237.6100
discoverysciences.AU@
grace.com

GRACE®, GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. REVELERIS® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. ALLTIMA™ and REVEALX™ are trademarks of Alltech Associates, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. © Copyright 2011 W. R. Grace & Co.-Conn. All rights reserved.