



Increase Sample Purity with GraceResolv™ High Resolution Flash Cartridges

**Romulus Gaita
Paul Garms
Scott Anderson**

Grace Davison Discovery Sciences
2051 Waukegan Rd.
Deerfield IL 60015 U.S.A.
Phone: 1-800-255-8324
Website: www.discoverysciences.com

GRACE

Abstract

The demand for fast and easy purification techniques is increasing. Flash chromatography can be a fast, yet crude methodology for preparative separations. To satisfy the need for flash cartridges with better resolution and higher loading capacity, GraceResolv™ high-resolution flash cartridges have been developed by Grace Davison Discovery Sciences. The silica packed in GraceResolv™ cartridges was specifically designed (patent pending) to increase cartridge efficiency (plates/ meter) and reduce unwanted mixed-mode interactions thereby increasing resolution of closely eluting compounds. This results in an increase in sample recovered at the desired purity level. The increase in resolution also enables a significant increase in loading capacity, over two times greater in some cases, versus traditional flash cartridges. By purifying greater amounts on a given cartridge size, users can reduce purification costs and increase productivity.




Introduction

Flash chromatography, also known as medium pressure chromatography was popularized several years ago as an alternative to slow and often inefficient gravity-fed chromatography. Today, flash chromatography has advanced considerably by introducing automated instruments using pre-packed flash chromatography cartridges eliminating the time-spent packing. The one area of flash chromatography that has seen little advancement is the heart of the separation, the silica.

Most cartridges today are packed with the same 40-60µm irregular shaped silica that was used for self-packed glass columns in the 1970s. The benefit of this type of silica is that it is inexpensive and generates low backpressure however it only provides a marginal level of resolution. Grace, as a silica manufacturer, and with extended chromatography knowledge is applying new technologies to the field of flash chromatography.

The silica used in GraceResolv™ cartridges is smaller in particle size, centered around 40um,with a more narrower particle size distribution than standard flash silica. The smaller particle size increases the efficiency of the cartridge resulting in higher resolution. While the narrow particle size distribution (**Figure A**) reduces fines (small particles less than 5um) which keeps backpressure low. The GraceResolv™ silica also undergoes a proprietary process to remove metals from the silica surface resulting in less peak tailing especially for metal chelator compounds (**Figure B**).

Figure A – GraceResolv™ Particle Size



A smaller particle size increases resolution. Tighter particle size distribution and minimal fines (small particles) keep backpressure low. Measured by Malvern® Mastersizer® 2000.

Figure B – GraceResolv™ Silica Purity

Percentage of Metal Content				
Silica	Fe ₂ O ₃	Na ₂ O	Al ₂ O ₃	SO ₄
GraceResolv™	0.001	0.01	0.02	0.001
Competitive Standard	0.005	0.07	0.05	0.01

Grace's flash silica undergoes a proprietary process to reduce metals from the silica surface by over 50%. Metals cause peak tailing and reduce resolution especially for metal chelator compounds.

Purpose


To demonstrate the advantages of higher resolution flash cartridges we compared GraceResolv™ cartridges to leading competitive cartridges. Experiments designed to test sample resolution, loading capacity, recovery, and throughput were performed.

Results and Discussion

Figure 1 demonstrates GraceResolv™ cartridges superior resolution of closely eluting impurities, separating two additional impurity compounds from the target drug and better resolving the main impurity verses leading competitor cartridges. This can dramatically increase the yield of high purity fractions as demonstrated in **Figure 2**, the purification of α from δ tocopherol. The increase in resolution also allows for higher loading capacity while maintaining base-line separation. This can enable chemists to use smaller cartridge sizes to purify the same amount of material typically purified on larger competitor cartridges reducing consumable costs (**Figure 3**). The combination of higher resolution and loading capacity also can enable chemists to increase flow rate and still maintaining baseline resolution, saving time (**Figure 4**).

Figure 1

Stage 1 Anti-obesity compound purified from synthetic byproducts (Data courtesy of Cambridge Biotechnology Ltd, part of the Biovitrum Group)



Increase sample purity using GraceResolv™ cartridges

Comparison of GraceResolv™ 40g flash cartridges to Competitor A 40g cartridge and Competitor B 50g cartridge demonstrates the increase in resolution possible with GraceResolv™ cartridges. The GraceResolv™ cartridge was able to separate 3 distinct impurity peaks and trigger 4 different fraction sets. The Competitor A cartridge could only separate 1 distinct impurity peak and trigger 2 fraction sets, while the Competitor B cartridge could not trigger a separate fraction collection from the main peak.

Purification of α -tocopherol from δ -tocopherol





Figure 3

Loading comparison of a 3.35g Phthalate mixture: GraceResolv™ 40g cartridge vs Competitor 80g cartridge.




Reduce costs using smaller GraceResolv™ cartridges

3.36g (8.4%) loading on 40 gram GraceResolv™ cartridge has a better resolution then an identical loading of 3.36g (4.2%) on an 80 gram Competitor A cartridge. This indicates a >2X loading capacity and the ability to down size the cartridge reducing cartridge costs and solvent consumption.

Figure 4

Flow rate comparison in the purification of α from δ -tocopherol.



Save time with faster purifications using GraceResolv™ cartridges

Higher resolution of GraceResolv™ cartridges enables the flow rate to be increased 2X while maintaining base-line separation. For this application, that results in over a 40% time savings per run.

Conclusion

GraceResolv™ cartridges are packed with a new silica (patent pending) that results in superior resolution and loading capacity versus the leading competitors. Increased resolution and loading capacity give GraceResolv™ cartridges:

- 3X increase in pure compound recovered
- 2X faster purification maintaining baseline resolution
- 30% reduction in consumable costs

ALLTECH® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. BRINGING SEPARATIONS TOGETHER™ is a trademark of Alltech Associates, Inc. GRACERESOLV™ and GRACE DAVISON DISCOVERY SCIENCES™ are trademarks of W. R. Grace & Co.-Conn. GRACE® and GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. COMBIFLASH® is a registered trademark of Teledyne Isco, Inc. MALVERN® and MASTERSIZER® are registered trademarks of Malvern Instruments, Limited. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership.

Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. © Copyright 2008 Alltech Associates, Inc. All rights reserved.

The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. Grace reserves the right to change prices and/or specifications without prior notification. Printed in the USA.