

# How to Choose a GC Capillary Column...

When selecting a capillary column for an application, there are four basic parameters that need to be considered:

- |          |                          |
|----------|--------------------------|
| <b>1</b> | Stationary Phase         |
| <b>2</b> | Column Internal Diameter |
| <b>3</b> | Film Thickness           |
| <b>4</b> | Column Length            |


## Choosing a Stationary Phase

Select a relatively non-polar stationary phase. SGE recommends a 5% phenyl phase which can be used for over 80% of analyses. It offers just enough polarity to be useful, but will still essentially separate components on boiling point.

Use a specialized phase when recommended by standard methods, for example BPX70 for Fatty Acid Methyl Ester (FAME) analysis.

**Figure 1** demonstrates the effect of different phase polarities on the same analysis. BP1 is 100% dimethylpolysiloxane, which is non-polar, while BP20 is a moderately polar polyethylene glycol phase.

**Figure 1. Effect of Phase Type**



## Stationary Phases:

- **BP** Industry Standard Phases
- **BPX** MS Premium Phases
- **HT** Very high temperature Phases




## Choosing the Column Internal Diameter

Of the four parameters to be considered when selecting a capillary column, internal diameter is often the easiest to select.

- A 0.25mm ID column is a good general selection for most applications.
- Columns with a larger bore (e.g. 0.53mm ID) will be more rugged and have a larger sample capacity.

**Figure 2** shows the effect of changing the internal diameter of a column will have on an analysis. Decreasing the diameter will increase resolution(R), or you will achieve the same resolution with a FASTER run time.

**Figure 2. Effect of Column Diameter**




## Film Thickness

Increasing the film thickness will increase retention and improve resolution (figure 3). Thicker films will reduce the possibility of large overloaded peaks co-eluting with other compounds.

A large range in solute concentrations will require a thicker film, while if the separation of solutes is sufficient and co-elution is unlikely, a thinner film can be used. Thicker films should be used for analysis of volatiles.

A change in the film thickness of a column will directly affect the temperature at which a compound will elute, the thicker the film, the greater the retention and thus the higher the elution temperature.

**Figure 3** shows the effect varying a column's film thickness will have on an analysis.





## Column Length

Length is often the last parameter to be considered when selecting a column.

- The longer the column; the greater the efficiency.
- Select the shortest column length which will provide the necessary resolution.
- If the maximum available column length gives inadequate resolution, consider a change in stationary phase or phase ratio.
- Resolution is proportional to the square root of column efficiency, doubling the column length will only increase the resolving power of a column by approximately 40%.

**Figure 4** shows how various column lengths affect an analysis.



**EXPERT Tips**

- 1 If peak shape deteriorates, replace your inlet liner and cut approximately 30cm from the front end of the capillary column.
- 2 When removing a column from the GC, seal the ends with an old septum to preserve column for future use.

*For further information, contact your nearest SGE office or distributor.*