

Columns for Silica-based Reversed Phase Chromatography (Other Columns)

Please refer to “Comparison of the Features Among Shodex Reverse Phase Chromatography (RPC) Columns” on page 6 and 7 for features.

Standard columns

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (μm)	Carbon Load (%)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F6650052	Silica 5C8 4D	≥ 9,000	Octyl	5	10	100	4.6 × 150	H ₂ O/CH ₃ OH=34/66
F6650053	Silica 5C8 4E	≥ 15,000	Octyl	5	10	100	4.6 × 250	H ₂ O/CH ₃ OH=34/66
F6650054	Silica 5C4 4D	≥ 9,000	Butyl	5	7	100	4.6 × 150	H ₂ O/CH ₃ OH=45/55
F6650055	Silica 5C4 4E	≥ 15,000	Butyl	5	7	100	4.6 × 250	H ₂ O/CH ₃ OH=45/55
F6650058	Silica 5CN 4D	≥ 7,000	Cyanopropyl	5	–	100	4.6 × 150	H ₂ O/CH ₃ OH=60/40
F6650059	Silica 5CN 4E	≥ 12,000	Cyanopropyl	5	–	100	4.6 × 250	H ₂ O/CH ₃ OH=60/40
F6650062	Silica 5NPE 4D	≥ 8,000	Nitrophenylethyl	5	–	100	4.6 × 150	H ₂ O/CH ₃ OH=45/55
F6650063	Silica 5PYE 4D	≥ 7,000	Pyrenylethyl	5	–	100	4.6 × 150	H ₂ O/CH ₃ OH=30/70

Base Material : Silica

Preparative columns * Preparative columns are made to order.

Product Code	Product Name	Plate Number (TP/column)	Particle Size (μm)	Column Size (mm) I.D. x Length	Standard Column
F7560062	Silica 5C8 10E	≥ 15,000	5	10.0 × 250	5C8
F7560063	Silica 5C8 20E	≥ 15,000	5	20.0 × 250	5C8
F7560054	Silica 5C4 10E	≥ 15,000	5	10.0 × 250	5C4
F7560055	Silica 5C4 20E	≥ 15,000	5	20.0 × 250	5C4

Columns for Silica-based HILIC and Normal Phase Chromatography

Features

- 5SIL**
- Uses high purity silica (99.99% or higher)
 - Suitable for normal phase analysis using a nonpolar organic solvent
 - Corresponds to USP L3

- 5NH**
- Suitable for saccharides analysis by hydrophilic interaction chromatography (HILIC)
 - Corresponds to USP L8

Standard columns

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (μm)	Carbon Load (%)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F6650050	Silica 5SIL 4D	≥ 9,000	–	5	–	100	4.6 × 150	C ₆ H ₁₄ /C ₂ H ₅ OH=95/5
F6650051	Silica 5SIL 4E	≥ 15,000	–	5	–	100	4.6 × 250	C ₆ H ₁₄ /C ₂ H ₅ OH=95/5
F6650060	Silica 5NH 4D	≥ 5,000	Aminopropyl	5	–	100	4.6 × 150	H ₂ O/CH ₃ CN=5/95
F6650061	Silica 5NH 4E	≥ 8,000	Aminopropyl	5	–	100	4.6 × 250	H ₂ O/CH ₃ CN=5/95


Base Material : Silica

Preparative columns * Preparative columns are made to order.

Product Code	Product Name	Plate Number (TP/column)	Particle Size (μm)	Column Size (mm) I.D. x Length	Standard Column
F7560050	Silica 5SIL 10E	≥ 15,000	5	10.0 × 250	5SIL
F7560051	Silica 5SIL 20E	≥ 15,000	5	20.0 × 250	5SIL
F7560060	Silica 5NH 10E	≥ 8,000	5	10.0 × 250	5NH
F7560061	Silica 5NH 20E	≥ 8,000	5	20.0 × 250	5NH

Aflatoxins


TFA method



Sample : 5µg/L each, 20µL
 1. Derived Aflatoxin G₁
 2. Derived Aflatoxin B₁
 3. Aflatoxin G₂
 4. Aflatoxin B₂

Reference :
 Shoku-An No. 0816-1
 (August 16, 2011, Japan)
 "Test Methods Related to
 Total Aflatoxin" in Notice

Column : Shodex Silica C18M 4E
 Eluent : H₂O/CH₃CN/CH₃OH=60/10/30
 Flow rate : 1.0mL/min
 Detector : Fluorescence (Ex. : 365nm, Em. : 450nm)
 Column temp. : 40°C




Sample : 50µL
 1. Aflatoxin M₁ 1µg/L

Reference :
 Shoku-An No. 0723-5
 (July 23, 2015, Japan)

Column : Shodex Silica C18M 4E
 Eluent : H₂O/CH₃CN=75/25
 Flow rate : 1.0mL/min
 Detector : Fluorescence (Ex. : 365nm, Em. : 435nm)
 Column temp. : 40°C

Trichothecene mycotoxins




Sample : 20µL
 1. Nivalenol (NIV)
 2. Deoxynivalenol (DON)
 * 1mg/L each of DON and
 NIV were added to wheat sample

Reference :
 Shoku-An No. 0717001
 (July 17, 2003, Japan)
 Test method of deoxynivalenol

Column : Shodex Silica C18M 4E
 Eluent : H₂O/CH₃CN/CH₃OH=90/5/5
 Flow rate : 1.0mL/min
 Detector : UV (220nm)
 Column temp. : 40°C


Ochratoxin



Sample : 10µL
 1. Ochratoxin A 50µg/L

Column : Shodex Silica C18M 4D
 Eluent : H₂O/CH₃COOH/CH₃CN=43/2/55
 Flow rate : 1.0mL/min
 Detector : Fluorescence (Ex. : 333nm, Em. : 460nm)
 Column temp. : 40°C


Gingerol and shogaol



Sample : 0.1mg/mL each, 10µL
 1. 6-Gingerol
 2. 6-Shogaol

Column : Shodex Silica C18M 4D
 Eluent : (A) ; H₂O (B) ; CH₃CN
 Linear gradient : (B%) 40% to 70% (15min)
 Flow rate : 1.0mL/min
 Detector : UV (280nm)
 Column temp. : 40°C

Chlorogenic acid




Sample : 0.1mg/mL each, 10µL
 1. Chlorogenic acid
 2. Caffeine
 3. Caffeic acid

Column : Shodex Silica C18M 4D
 Eluent : 20mM H₃PO₄ aq. /CH₃OH=70/30
 Flow rate : 1.0mL/min
 Detector : UV (280nm)
 Column temp. : 30°C

Catechins


Sample : 10 μ g/mL each, 10 μ L
 1. Epigallocatechin
 2. Catechin
 3. Epigallocatechingallate
 4. Epicatechin
 5. Epicatechingallate
 6. Catechingallate



Column : Shodex Silica C18P 4D
 Eluent : (A) ; 20mM H₃PO₄ aq. (B) ; CH₃CN
 Linear gradient:
 (B%) 20% (0 to 5min), 20 to 40% (5 to 15min),
 40% (15 to 20min)
 Flow rate : 1.0mL/min
 Detector : UV (280nm)
 Column temp. : 30°C

Estradiols


Sample : 0.1mg/mL each, 10 μ L
 1. Estriol
 2. 17 β -Estradiol
 3. Estrone



Column : Shodex Silica C18P 4D
 Eluent : H₂O/CH₃CN=65/35
 Flow rate : 1.0mL/min
 Detector : UV (280nm)
 Column temp. : 30°C

Coenzyme Q10


Sample : 5 μ L
 1. Coenzyme Q10 1mg/mL



Column : Shodex Silica C18P 4D
 Eluent : CH₃OH/C₂H₅OH=13/7
 Flow rate : 1.2mL/min
 Detector : UV (275nm)
 Column temp. : 35°C

Benzylpyridine isomers


Sample :
 1. 2-Benzylpyridine
 2. 3-Benzylpyridine
 3. 4-Benzylpyridine



Column : Shodex Silica 5PYE 4D
 Eluent : 20mM KH₂PO₄ aq./CH₃OH=40/60
 Flow rate : 1.0mL/min
 Detector : UV (254nm)
 Column temp. : 30°C

Dinitronaphthalene isomers


Sample :
 1. Naphthalene
 2. 1,5-Dimethylnaphthalene
 3. 1,5-Dinitronaphthalene
 4. 1,8-Dinitronaphthalene



Column : Shodex Silica 5NPE 4D
 Eluent : H₂O/CH₃OH=30/70
 Flow rate : 1.0mL/min
 Detector : UV (254nm)
 Column temp. : 30°C

Simultaneous analysis of vitamin E homologs

Sample : 20 μ L
 1. α -Tocopherol 5 μ g/mL
 2. α -Tocotrienol 10 μ g/mL
 3. β -Tocopherol 5 μ g/mL
 4. γ -Tocopherol 5 μ g/mL
 5. γ -Tocotrienol 10 μ g/mL
 6. δ -Tocopherol 5 μ g/mL
 7. δ -Tocotrienol 10 μ g/mL



Column : Shodex Silica 5SIL 4D
 Eluent : n-Hexane/Isopropanol/Acetic acid=1000/6/5
 Flow rate : 1.0mL/min
 Detector : Fluorescence (Ex. : 298nm, Em. : 325nm)
 Column temp. : 30°C