

Inserts

Our recommended vials, caps and inserts are specifically designed and made for the autosamplers listed below.

Please see our dedicated autosampler pages for more information.

Ordering Information

PART NO. (100pk)	DESCRIPTION	FITS AUTOSAMPLER VIAL
V0300P-1232	300µL Conical, Ultraclean Polypropylene, Snap/Crimp Top <i>Ultra-Clean For Bio Assays</i>	
I025BS-629	250µL Glass Conical Insert with Bottom Spring	Agilent™ <ul style="list-style-type: none"> • 2mL Crimp, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 9mm thread • 2mL Snap, 12 x 32mm, 11mm opening Perkin Elmer Bruker Waters Alliance 2690 Waters Wisp <ul style="list-style-type: none"> • 2mL Snap, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 10mm thread • 2mL Screw, 12 x 32mm, 10mm thread • 1mL Shell, 8 x 40mm
I025-630	350µL Glass Flat Bottom Insert	Agilent™ <ul style="list-style-type: none"> • 2mL Crimp, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 9mm thread • 2mL Snap, 12 x 32mm, 11mm opening Perkin Elmer Bruker Waters Alliance 2690 <ul style="list-style-type: none"> • 2mL Snap, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 10mm thread • 2mL Screw, 12 x 32mm, 10mm thread
I025CP-631	250µL Conical Point Glass Insert, 6 x 31mm Stand alone Insert	Agilent™ <ul style="list-style-type: none"> • 2mL WM Crimp Vials
I025GF-631	250µL Glass BM Hanging Conical Point Insert	Agilent™ <ul style="list-style-type: none"> • 2mL Crimp, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 9mm thread • 2mL Snap, 12 x 32mm, 11mm opening Perkin Elmer Bruker Waters Alliance 2690 <ul style="list-style-type: none"> • 2mL Snap, 12 x 32mm, 11mm opening • 2mL Screw, 12 x 32mm, 10mm thread • 2mL Screw, 12 x 32mm, 10mm thread
I01-531	100µL Glass LV Insert, Conical Point, No Spring	Beckman/Hitachi Perkin Elmer Bruker <ul style="list-style-type: none"> • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread
I01BS-530 I01PBS-530	100µL Glass LV Insert with Bottom Spring 100µL PP LV Insert with Bottom Spring	Beckman/Hitachi Perkin Elmer Bruker Waters Wisp <ul style="list-style-type: none"> • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread • 1mL Shell, 8 x 10mm
I025-531	250µL Glass Flat Bottom Insert	Beckman/Hitachi Perkin Elmer Bruker <ul style="list-style-type: none"> • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread • 2mL Screw, 12 x 32mm, 8mm thread
I03-638X	250µL Glass Conical Insert, Requires Metal Spring	Waters Wisp <ul style="list-style-type: none"> • 4mL Screw, 15 x 45mm, 13mm thread
I03-639	250µL Glass Mandrel Point Insert, Requires Metal Spring	Waters Wisp <ul style="list-style-type: none"> • 4mL Screw, 15 x 45mm, 13mm thread
I025-629	250µL Mandrel Point Insert, Requires Metal Spring	Waters Wisp <ul style="list-style-type: none"> • 1mL Shell, 8 x 40mm
I02BS-536	250µL Glass Conical Point Insert with Bottom Spring	Waters Wisp <ul style="list-style-type: none"> • 1mL Shell, 8 x 40mm

