

EXtrelut® Columns

Liquid-liquid extraction in its most effective form



EXtrelut® NT sorbent columns simplify liquid-liquid extraction by replacing separation funnels. Using a single step is more efficient and saves solvent, material and time in contrast to classical funnel separation. With its intuitive working principle, higher recovery and cleaner extraction can be achieved. The aqueous sample is simply applied to the EXtrelut® NT sorbent, which distributes itself as a thin film over the chemically inert matrix and thus acts as a stationary phase.

Features & Benefits

- Save solvents
- Easy to use
- Highly efficient

Applications

Liquid-Liquid Extraction

Ordering Information

Description	Qty/Pk	Catalogue No.
EXtrelut® NT1	100 columns	1.15094.0001
EXtrelut® NT3	50 columns	1.15095.0001
EXtrelut® NT20	25 columns	1.15096.0001
EXtrelut® NT refill packs	50 bags	1.15093.0001

For more information visit: www.merckmillipore.com/extrelut

Filter Forceps

For safe handling of filter membranes



Ensure your filtration membranes are not damaged during handling by using these blunt-ended forceps.

Features & Benefits

- Beveled, unserrated tips to protect delicate membrane surfaces
- May be sterilized by autoclaving or flame-sterilization

Applications

Gentle Handling of Membranes

Ordering Information

Description	Qty/Pk.	Catalogue No.
Filter Forceps, blunt end, stainless steel	3	XX6200006P

For more information visit: www.merckmillipore.com/LabHardware