


InertCap™ for Amines

- Cross-Linked
- Excellent for amine compounds from C2 ~ C10 analyses
- Ideal for the simultaneous analyses of mixed sample such as alcohol etc.
- Equivalent : CP-Volamine

InertCap for Amines shows excellent performances for the amine compounds from C2 ~ C10. Basic compounds can be perfectly eluted without being adsorbed in the column. Unlike other amines special columns, InertCap for Amines can simultaneously elute general solvents such as alcohol due to our unique neutral inert treatment technique.

■ Solvent and amine mixture analyses



【InertCap™ for Amines】

ID (mm)	Length (m)	Thickness (μm)	Max. Temperature (°C)	Cat.No.
0.32	15	—	iso.265-prog.300	1010-69229
	30	—	iso.265-prog.300	1010-69249
	60	—	iso.265-prog.300	1010-69269