



SPECIALTY CARTRIDGES

FOR ENVIRONMENTAL EXTRACTION AND CLEAN UP



8

www.unitedchem.com

PRODUCT NAME	DESCRIPTION	PART NUMBER	AMT. SORBENT/ TUBE VOL.	UNITS
CLEAN-ELUTE™	A diatomaceous earth matrix capable of being used within a pH range of 1-13.	CLEAN-ELUTE	25,000mg/200mL	108
Florisil® PR	Pesticide Residue grade Florisil®	EUFLS1M6	1000mg/6mL	30
Florisil® A	Equivalent to PR grade, but with a smaller particle size. Preferred by many environmental testing labs.	EUFLSA1M6	1000mg/6mL	30
ENVIRO-CLEAN® EPH Fractionation	(Gravity Flow) Developed for fractionation of MA EPH.	XRSIHT13M15 CUSILHT15M25	3000mg/15mL 5000mg/25mL	24 20
ENVIRO-CLEAN® C18	Optimized for non polar analytes, including pesticides, PCBs, PAHs, and formaldehyde.	EEC181M6	1000mg/6mL	30
ENVIRO-CLEAN® C18 Polar	For non polar and moderately polar analytes.	EUC181M6	1000mg/6mL	30
ENVIRO-CLEAN® 548	For EPA Method 548.1 Endothall extraction.	EC548006	6mL	30
ENVIRO-CLEAN® PS-DVB	EPA Method 526, 528, 529 and 8330	ECDVB156	500mg/6mL	30
ENVIRO-CLEAN® 521 and 522	Activated carbon for EPA Methods 521	EU52112M6	2000mg/6mL	30
ENVIRO-CLEAN® 523	Graphitized Carbon for Method 523	EC5232506	250mg/6mL	30
ENVIRO-CLEAN® 525	Specially designed C18 for EPA Method 525.2	EC525006	1500mg/6mL	30
ENVIRO-CLEAN® (GCB) 90² m/g Surface Area	For applications requiring graphitized carbon black with lower surface area.	EC535156	500mg/6mL	30
ENVIRO-CLEAN® 537	For EPA Method 537 Perfluorinated Alkyl Acids	ECDVB156P	500mg/6mL	30

Florisil® is a registered trademark of U.S. Silica.



UCT, Inc. | 2731 Bartram Rd | Bristol, PA 19007 | USA
P. 800.385.3153 | 215.781.9255 | F. 215.785.1226
www.unitedchem.com



ENVIRO



SPECIALTY CARTRIDGES AND DISKS

FOR ENVIRONMENTAL EXTRACTION AND CLEAN UP



6 mL



PRODUCT NAME	DESCRIPTION	PART NUMBER	AMT. SORBENT/ TUBE VOL.	UNITS
ENVIRO-CLEAN® Silica Gel Cartridge	For silica gel clean-up applications	EUSILMSSM26	1000mg silica/200mg muffled sodium sulfate anhydrous/6 mL	30
ENVIRO-CLEAN® PAH Extraction Salts	For AOAC Method Extraction of PAHs from Gulf of Mexico Seafood	ECQUUS2-MP	4000mg magesium sul- fate anhydrous/2000mg sodium chloride	50
ENVIRO-CLEAN® PAH	Centrifuge Tubes for PAH extraction	ECPAHR50CT	50mL	50
ZERO-BLANK™ FILTER	To prevent airborne SPE contamination	ECBLANK	83mL	6
ENVIRO-CLEAN® C08	C08 for EPA Method 549 Diquat and Paraquat	EEC08156	500mg/6mL	50
SODIUM SULFATE ANHYDROUS Drying Cartridge	A convenient product for the drying of extracts prior to concentration and analysis.	ECSS15M6	5000mg/6mL	30
HIGH LIPID CLEAN-UP	For the clean-up of agricultural matrices containing high levels of lipids. See application note.	ECPSAC1856	500mg C18 on top / 500mg PSA on bottom	30
PS-DVB Push-Thru Cartridge	DVB in a Push-Thru cartridge	ECDVBMC	160mg	25
Sodium Sulfate Push Thru Cartridge	A convenient product for the drying of extracts prior to concentration and analysis	ECSSLC	3000mg	100
Endcapped C18 Push-Thru Cartridge	Endcapped C18 in a Push-Thru cartridge	EEC1815MC-P	500mg	50

