

XtrackT[®] GRAVITY FLOW SPE COLUMNS

XtrackT[®] large particle bonded phases allow for uniform gravity flow for most blood and urine samples. A single column provides extraction for a broad spectrum of compounds with selective elution of acid neutrals, steroids and bases. XtrackT[®] large particle (80-200 µm) silica gels are available with hydrophobic, hydrophilic, ion exchange or copolymeric phases, including DAU mixed mode. XtrackT[®] is recommended for viscous sample matrices or for gravity flow applications.



XtrackT[®]

GRAVITY FLOW COLUMNS

GRAVITY FLOW XtrackT[®] DAU SORBENT (XRDAH)

Organic Loading = 12.4%
Surface Area = 500 m²/g

Average Pore Size = 60 Å
Pore Volume = 0.77 cm³/g

Cation Exchange = 0.077 meq/g

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
3	150	50	No	XRDAH(150)03
3	200	50	No	XRDAH203
3	200	50	Yes	XCDAH203
3	300	50	No	XRDAH303
3	500	50	No	XRDAH503
3	500	50	YES	XCDAH503
6	200	50	NO	XRDAH206
6	500	50	NO	XRDAH506
10	130	50	NO	XRDAH13Z
10	130	500	NO	XRDAH13Z-D
10	200	50	NO	XRDAH20Z
10	200	50	YES	XCDAH20Z
10	500	50	NO	XRDAH50Z
15	500	50	NO	XRDAH515
15	500	50	YES	XCDAH515
15	1000	50	NO	XRDAH15

18

GRAVITY FLOW XtrackT[®] ENDCAPPED C18 COLUMNS (XRODH)

Organic Loading = 21.5%
Surface Area = 500 m²/g

Average Pore Size = 60 Å
Pore Volume = 0.77 cm³/g

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
3	500	50	NO	XRODH503
3	500	500	NO	XRODH503-D
3	500	50	YES	XCODH503
6	500	50	NO	XRODH506
6	500	50	YES	XCODH506
6	1000	50	NO	XRODHM06
15	500	50	NO	XRODH515
15	1000	50	NO	XRODHM15
25	5000	20	NO	XRODH5M25
75	10000	10	NO	XRODH10M75

GRAVITY FLOW XtrackT[®] BENZENESULFONIC ACID SORBENT (XRBSH)

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
10	500	50	NO	XRBSH50Z
15	500	50	NO	XRBSH515

GRAVITY FLOW XtrackT[®] CARBOXYLIC ACID SORBENT (XRCCH)

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
3	500	50	NO	XRCCH503
6	500	50	NO	XRCCH506
15	1000	50	NO	XRCCHM15

GRAVITY FLOW XtrackT[®] PROPYLSULFONIC ACID SORBENT (XRPCH)

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
3	200	50	NO	XRPCH203
3	500	50	NO	XRPCH503
10	500	50	NO	XRPCH50Z

GRAVITY FLOW XtrackT[®] PRIMARY/SECONDARY AMINE SORBENT (XRPSH)

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
3	300	50	NO	XRPSH303

GRAVITY FLOW XtrackT[®] HEAT TREATED SILICA SORBENT (XRSIHT)

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	CLEAN-THRU [®] Tips Provided	Part Number
10	500	50	NO	XRSIHT50Z
15	3000	24	NO	XRSIHT13M15

*XRSIHT13M15 also comes with Flange Caps and Luer Tips