



オレイン酸とエライジン酸の分離

Cis-trans isomers of 9-cotadecenoic acid

SunShell C18 2.6 μm , 150 x 4.6 mm i.d.



Column: SunShell C18, 2.6 μm 150 x 4.6 mm
Mobile phase: Acetonitrile:water:formic acid=90:10:0.05
Flow rate: 1.0 mL / min
Temperature: 40 $^{\circ}\text{C}$
Detection: UV@215 nm
Sample: 1 = Oleic acid, 2 = Elaidic acid